 **GB** B. Braun Medical Ltd.
Sheffield, S35 2PW
United Kingdom

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





226306*
226308*
226310*
226312*
226314*
226316*

* For each country a specific letter is added at the end of the reference

Actreen® Glyc Set

EN Instructions for use

Carefully read the Instructions for use before usage

Device description

Actreen® Glys Set is a pre-lubricated, sterile and ready to use catheter connected to a graduated bag. The bag includes a non-return valve.

Composition

Actreen® Glys Set is made of PVC (catheter), PE (pouch), the lubricant is water and glycerin based and the primary packaging is made of PE/PET and PGL.

Not made with DEHP and Latex.

Sterility

Sterilized using Ethylene Oxide.

Intended purpose

The intended use (purpose) of Actreen® Glys Set is to provide intermittent urinary catheterisation of the bladder, by introducing the catheter through the urethra to drain the urine from the bladder.

Indication

Actreen® Glys Set is indicated for patients who need intermittent urinary catheterisation due to temporary or chronic urine retention or voiding dysfunction.

Patient population

Actreen® Glys Set can be used in male and female adult or paediatric patients with a body weight superior or equal to 3,5 kg.

The responsible healthcare professional shall prescribe the appropriate size, tip and length of the catheter as well as the frequency of intermittent catheterisation based on patient's age, gender, characteristics, and conditions (body profile in general, pregnancy, paediatric population) according to state-of-the-art medical knowledge.

For pregnant women, the responsible healthcare professional shall assess the appropriate length and diameter of the catheters, since in pregnancy the length of the urethra is altered.

For paediatric patients, the capability to perform intermittent self-catheterisation should be determined by the responsible healthcare professional according to the state-of-the-art medical knowledge.

Intended User

The intended users are healthcare professionals, caregivers, and patients.

Caution: The responsible healthcare professional must ensure that the patients and caregivers are duly trained and instructed in intermittent (self) catheterisation with the catheter, as use of the device by an untrained user may cause trauma or compromise the sterility of the product, potentially leading to infection.

Intended Use Environment

Home, daily life activities environments, hospitals and healthcare facilities.

Contraindications

Contraindication related to the product :

Known hypersensitivity reaction to any material used (see composition).

Absolute contraindications to intermittent catheterisation:

High intravesical pressure. Urinary fistula. Mitrofanoff procedure. Priapism. After urethral surgery. Suspected urethral tumor. Pelvic / urethral injury.

Relative contraindications to intermittent catheterisation with special cautions:

Please contact your healthcare professional before using this product if any of the following conditions are present: Abnormal urethral anatomy (e.g., stricture), false passage, bladder neck obstruction, diseases of penis (e.g., injury, tumor, or infections), after gender reassignment surgery, tendency for autonomic dysreflexia with bladder filling/urethral instrumentation (despite treatment), urethral bleeding, prostatic or bladder neck surgery, female genital mutilation, stent or artificial prosthesis in the urogenital tract, bladder capacity < 200 ml, high fluid intake regimen.

Contraindication for intermittent self-catheterisation:

Inability to self-catheterize in the absence of an appropriately trained caregiver or attendant for patient with poor manual dexterity and cognitive impairment.

Warning: If you have any of the above-named contraindications or are in doubt, please contact your responsible physician immediately before the next catheterisation on how to proceed.

The responsible physician must define contraindications for an intermittent catheterisation based on patient characteristics and conditions according to state-of-the-art medical knowledge.

Residual Risks

Residual risks are included in other chapters (e.g., warnings etc.).

Side effects

Side effects related to the product

Unknown hypersensitivity reaction to any material used leading to anaphylactic or anaphylactoid reactions.

Side effects related to intermittent catheterisation

Catheter-associated Urinary Tract Infection (CAUTI), in some instance associated with pyelonephritis. Urethritis / epididymo-orchitis / prostatitis. Trauma of the urinary tract including urethral trauma / urethral bleeding / false passage / urethral strictures (predominantly in men) / meatal erosion and stenosis / bladder perforation. Discomfort, soreness, pain. Bladder calculus, stone formation.

Warnings: If severe hypersensitivity reactions occur, remove the product immediately, discontinue its use and contact a physician. If you experience reduced or interrupted urine flow or that the catheter does not drain the bladder properly, contact a physician as it may indicate that you should use another size or type of catheter. If you have symptoms of CAUTI or infection of urogenital tract: fever, chills, pyuria, dysuria, urgency, cloudy urine, suprapubic pain, discomfort in emptying the bladder, blood in the urine, increased bladder spasms, autonomic dysreflexia, mucus, contact a healthcare professional.

If you have any signs of trauma, discomfort when emptying the bladder, bleeding, a frequent urge to urinate or blood in the urine, discontinue use immediately and contact a healthcare professional. Ungentle insertion and removal can cause trauma and aggravate discomfort and pain. If you have doubts about the catheter handling, contact a healthcare professional.

Avoid introduction of pubic hair. Bladder calculus or stone formation are related to the introduction of pubic hair and presence of mucus.

Cautions

Fluid intake should be adequate unless contraindicated. Patient should avoid excessive fluid intake that may cause bladder overdistension.

Constipation may cause pressure preventing drainage from catheter: Maintain regular bowel function.

For all intended users, a trained period is necessary to be autonomous on catheterisation.

The healthcare professional must ensure that the patients and caregivers are duly trained and instructed in intermittent (self) catheterisation with the catheter, as use of the device by an untrained user may cause trauma or compromise the sterility of the product, potentially leading to infection.

Caution specific for healthcare professionals: Apply with caution if the patient produces urine with many particles clearly distinguishable by the naked eye, as it may lead to transient urine retention.

Catheterisation should take place as quickly as possible after opening the package as it may promote contamination or compromise product characteristics (e.g., drying lubricant).

Do not fold the catheter as it may potentially promote trauma to the urinary tract.

Avoid contact between the urinary meatus and the connector or the sheath as it may potentially promote trauma to the urinary tract.

Urine splash may occur in case of intense pressure on the bag or bag carried upside-down.

Warning

Any contamination of the catheter needs to be avoided by proper handling.

Do not use a female size catheter for an adult male patient.

Do not use a male Tieman catheter for female patient.

Cleaning, disinfection, resterilisation may damage the product and creates a potential risk to the patient. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

Do not use the catheter if you observe any damage (e.g., kinked) to, blockage of, or foreign bodies in the device.

Do not use if side eyes are missing as the device may not perform as intended.

For urethral use only.

Catheters must be removed immediately after completion of urinary drainage.

In case of suspected catheters knotting, please contact your healthcare professional immediately.

Do not use if package is damaged or after expiry date, as sterility of the product may be compromised, potentially leading to infection.

The scale printed on the bag gives an approximative volume of the urine content. Do not use for diuresis purpose.

Notice: There are additional, specific warnings in other sections of this Instruction for Use. These specific warnings must be considered as well.

Operating Instructions

1. PREPARATION:

•Carefully wash hands before and after catheterisation and clean the area surrounding the urethral opening (1A/1B).

If a caregiver is assisting catheterisation, local guidelines must be applied.

Before opening, remove the blister, remove the protective sleeve which protects the sheath and unfold the bag. Check the product is intact (no packaging damage, no visual defects, kink on tube...). Make sure the lubricant is evenly spread over the catheter. If not, spread the lubricant by touching the catheter over the packaging to allow an even distribution.

Unevenly spread lubricant may potentially promote trauma of the urinary tract.

During opening, keep the product horizontally and remove protective film slowly to avoid lubricant loss.

2. PRODUCT OPENING:

For safety and hygienic reasons, we advise you to use the product in shortest possible time after opening.

Caution: Pay attention no to put the catheter tube in contact with hand or objects/surfaces.

Warning: Wrong manipulation of the catheter could lead to contamination and/or impairment of functional capability.

Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

A. Peel open the outer packaging and remove the catheter (2).

B. Remove the protective sleeves from around the catheter (3). You should now be left with just the sheath surrounding the catheter itself. Unfold the bag.

C. Unfold the top of the sheath and tear open on the perforated line (4).

D. The catheter is sterile and lubricated. Push the tip of the catheter out of the sheath a short distance (5). Do not touch the tip. You can hold the catheter through the sheath.

3. CATHETERISATION:

Precaution: Pay attention not to fold the pouch/bag to allow a correct flow of urine inside.

MALE

- With one hand hold the penis upwards to reduce curvature of the urethra (6B (1)).
- With the other hand carefully introduce the catheter into the urethra until it reaches the bladder. Point your penis downwards in case of resistance. Continue inserting the catheter until urine starts to flow (6B(2)).
- Tiemann/Coudé catheter : In order to correctly guide the catheter, ensure to maintain the upward orientation of the tip all along the catheterisation.

Warning: The curved tip of the catheter shall point upwards. Wrong orientation of the tip may lead to trauma of the urinary tract.

FEMALE

With one hand spread the outer lips of the vulva to locate the urethral opening. With the other hand carefully introduce the catheter into the urethra until it reaches the bladder and the urine drains (6A).

MALE / FEMALE

- The urine drains into the bag (7A) (7B). Hold the catheter in place until the bladder has completely been drained. An eyelet is incorporated at the top of the bag, in order to hold it or hang it.
- Manual pressure on the lower abdomen towards the end of catheterisation will ensure that the bladder has been completely emptied.
- Immediately after the end of bladder emptying, withdraw the catheter.
- To stop the catheter from dripping, turn up the connector end of the catheter when removing.

The volume of urine evacuated could be identified-estimated on the graduated urine bag.

4. EMPTYING:

The tear-off upper corner of the bag may be removed (8).

The bag may then be emptied (9).

At this stage, a urine sample may be collected for analysis.

Duration of use

For single use only.

Warning: Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness, or death of the patient.

Catheter must be removed immediately after completion of urinary drainage.

Disposal

Discard the catheter according to the current local regulations (10).

Storage and handling conditions

Keep away from sunlight and heat. Keep dry. The product must be stored flat.

Warning: Do not store the product (unit) in contact with sharp materials.

Warning : the expiry date only applies when Actreen® Glys Set is correctly stored in an undamaged packaging.

Do not use after the expiry date.

Clinical benefits & Performances

Actreen® Glys Set allows the drainage of urine from the bladder through the urethra.

Notice to the user

If, during the use of this product or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

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