




AESCULAP® Safety Scalpels



Safe and sharp

AESCULAP® Safety Scalpels

Safe and sharp*



Features of the AESCULAP® Safety Scalpel in compliance with TRBA 250



"click" palpable

- The safety mechanism is part of the system and compatible with other accessories.
- Single-handed activation must be possible.
- ✓ The mechanism is activated automatically by sliding back the blade with one hand.
- The safety mechanism prevents reuse.
- ✓ The blade is safely locked in the disposal position.
- Using the safety product does not require a different technique of use.
- ✓ It is used in the same way as a standard single use scalpel.
- The activation of the safety mechanism must be clearly indicated (palpable or audible).
- ✓ An audible „click“ indicates that the safety locking mechanism has been activated, and this can be seen and felt immediately.

AESCULAP® Safety Scalpels

Safe and sharp*

Ordering Information



BA810SU
Fig. 10



BA811SU
Fig. 11



BA812SU
Fig. 12



BA815SU
Fig. 15



BA815-CSU
Fig. 15-C



BA818SU
Fig. 18



BA819SU
Fig. 19

Sales unit PAK = Package of 10 pieces



BA820SU
Fig. 20

BA821SU
Fig. 21

BA822SU
Fig. 22

BA823SU
Fig. 23

BA824SU
Fig. 24

BA825SU
Fig. 25

BA836SU
Fig. 36

1/1

AESCULAP[®] – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 7461 95-0 | Fax +49 7461 95-2600 | www.aesculap.com

The main product trademark "Aesculap" is a registered trademarks of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.